



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SECA 207 Infant growth meter (Ex-Display)

- New, ex-display,
- Real photos of the product,
- Made in Germany,
- Seca 207 is an infant growth meter with extended caliper for precise measurement of a lying infant,
- It can be permanently bolted to a changing table, mounted on a wall or used as portable device,
- Large head and foot stops allow easy measurement of length,
- Folding zippers for easy storage
- Durability due to aluminum construction,
- Specification:
 - Measurement range: 0 - 99 cm,
 - Scale interval: 1 mm,
- Included:
 - Multilingual user manual (PDF),
 - Overall dimensions: 104.7 x 12 x 28.6 cm,
- Weight: 0.8 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

