



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



TELOS CPM Rehabilitation splint (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 45 W,
- TELOS CPM rail is used for the treatment of restrictions in the range of motion of the knee joint and also indirectly of the hip,
- Specifications:
 - Flexion angle adjustable in 4° increments from 0 - 98°,
 - Automatic length compensation, adaptation to the anatomical movement of the leg,
 - Can be used for both right and left leg,
 - Adjustable interval between flexion and extension in range: 1-99 seconds,
 - Adjustable speed in range: 60 - 100 %
 - Has memory - remembers last settings,
 - Emergency switch - the patient can stop the treatment at any time and start it again in the opposite direction,
- Control panel in German,
- Mobile on a cart,
- Includes:
 - User manual in German (PDF),
 - Power cord,
- Dimensions:
 - Rail: 66 x 30 x 35 cm,
 - Overall: 75 x 55 x 117 cm,
- Weight: 13 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

