



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ORTHOFEX MICRO-AIR Pedicure podiatric drill with dust extractor (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Netherlands,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- ORTHOFEX Micro-Air is a podiatric device with dust extraction, developed for podiatrists with professional pedicure in mind. Designed for continuous operation, allowing treating more patients. Adjustable speed with a maximum speed of 40,000 rpm. Suction ensures that the work surface is always clean, all dust is collected in the dust bag,
- Specifications:
 - Powerful motor - 40,000 rpm,
 - Variable speed control,
 - Right-left rotation,
 - Integral vacuum,
 - Low noise level,
 - Lightweight handle,
 - Rubber dust covers,
 - Replaceable dust bags,
 - No mill or drill included,
- Included:
 - Foot switch,
 - 2 x dust bags,
 - Additional fuse,
 - Power cord,
- Dimensions: 25 x 16 x 16 cm,
- Weight: 3.5 kg
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

