



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **DRAGER JULIAN Anesthesia machine (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in Germany,
- Power supply: 230 V,
- Emergency power supply: via 2 new 12 V 7 Ah batteries,
- Frequency: 50/60 Hz,
- Drager Julian is an inhalation anesthesia machine used in operating, preoperative and recovery rooms,
- Features two ports for vaporizers with common Dräger PLUG-IN connector systems, no vaporizers included - can be retrofitted at additional cost,
- Powered by gases: 2 x O<sub>2</sub>, 2 x N<sub>2</sub>O, Air,
- Bellows in Bottle type pneumatically driven respirator (electronically controlled),
- Ventilator measurement functions:
  - P<sub>MAX</sub> (V) - 10-70 mbar,
  - V<sub>T</sub> - 50-1400 ml,
  - f - 4-60 bpm,
  - Ti/Te - 1:4 to 2:1 with IPPV, 1:4 to 1:1 with PCV,
  - Tip/Ti - 0% - 50%,
  - PEEP - 0-15 mbar,
  - P<sub>insp</sub> (P) - 10-70 mbar,
  - Insp Flow - 10-75 L/min,
- Built-in LED monitor with O<sub>2</sub>, AIR, N<sub>2</sub>O gas measurement function, control panel and ventilator function display,
- Infinitely adjustable O<sub>2</sub>, AIR, N<sub>2</sub>O gas flow,
- Built-in suction unit powered from the central vacuum system,
- Airway pressure: from -10 to 99 mbar,
- Automatic recognition of anesthetic gases,
- Easy to use,
- Selectable menu in PL, EN, FR, JA, EL, HU, CA, TR, DE, IT, ES, NL, DA, NO, SV, PT,
- Keypad on the control panel in English,
- Built-in alarms, low and high pressure valves,
- Machine equipped with:
  - Expandable table with dimensions: 34 x 33.5 cm, max load 25 kg,
  - Drawer with dimensions: 30 x 44 x 18 cm,
  - Rail for attaching additional accessories,
  - Built-in suction unit powered from the central vacuum system,
  - Space for two gas cylinders,
  - Mobile on four casters including two with brakes,
- Includes:
  - Gas supply hoses with tips of customer's choice (AGA or DIN),
  - Brand new patient breathing system,
  - User manual in Polish (PDF),
- On request, at an additional cost, the possibility to:
  - Retrofit the device with a patient monitor/cardiac monitor from our offer,
  - Installing an anesthetic gas scavenging system (AGS, AGSS),
- Overall dimensions: 90 x 72 x 134 cm,
- Weight: 96 kg,
- In the domestic market (Poland), in the case of installation of the device within 30 days from the date of delivery of the equipment, installation and training are free of charge, after exceeding this time limit, these services will be charged,
- Warranty:
  - 6 months (for the domestic market - Poland),
  - 3 months for the international market,
- Financing option (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

