



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OMRA-WF Ophthalmic lens for laser treatments (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in the USA,
- OMRA-WF lens for trans-retinal photocoagulation in diabetic retinopathy. Excellent binocular resolution across the entire field of view. Allows the use of a large range of magnifications in slit lamps.
- Specifications:
 - Magnification: 0.68 x,
 - Magnification for laser beam: 1,50 x,
 - Contact area size: 15.5 mm,
 - Lens height: 28 mm,
 - Static field of view: 118°,
 - Dynamic field of view: 127°,
- Weight: 0.3 kg
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

