



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



AIR LIQUIDE GT35 Long-term storage tank for biological samples (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- Real photos of the product,
- Made in France, (2017),
- GT35 tank for long-term storage of biological samples in straws or vials features very low consumption of liquid nitrogen. Dedicated to storing a large number of samples in a small volume,
- Specification:
 - Made of aluminum,
 - Usable volume: 33.6 l,
 - Neck diameter: 50 mm,
 - Daily evaporation: 0.09 l,
 - Static holding time: 350 days,
 - Number of containers: 6,
 - Diameter of containers: 38 mm,
 - Height of containers: 280 mm,
 - Number of levels: 1,
 - Max capacity of 0.25 ml straws: 3120,
 - Max capacity of 0.5 ml straws: 1440,
 - Max capacity of 2 ml vials: 180,
- Has normal signs of use,
- Includes:
 - 6 280 mm containers,
- Dimensions:
 - Diameter: 47 cm,
 - Height: 66 cm,
- Weight:
 - Empty container: 15 kg,
 - Full tank: 43 kg,
- Valid technical passport issued,
- Warranty:
 - 6 months for domestic market - Poland,
 - 3 months for the international market,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

