



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MELAG 100 sealer (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Germany,
- Date of production: 2007,
- Power supply: 230 V,
- Frequency: 50/60 zHzHyHz,
- Power: 300 W,
- **Melag 100** is a paper and foil package sealing machine,
- Made of stainless steel,
- Immediately ready to work when turned on,
- An in-built cut-off knife,
- Smooth temperature control with a knob,
- Paper-foil sleeve holder,
- Pressure handle mounted on the right or left side,
- Seal length: 25 cm,
- Seal width: 1 cm,
- Overall dimensions: 35 x 40 x 33 cm,
- Weight: 5,5 kg,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (service report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

