



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEMMERT LABORATORY INCUBATOR BE 800 MODEL (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Germany
- Power supply: 230 V,
- Frequency: 50/60 zHzHyHz,
- Power: 2000 W,
- **Memmert BE 800** is a laboratory incubator with natural air circulation,
- It is used to sterilize medical materials by dry heating method with hot air under atmospheric pressure,
- PID electronic process controller with programmable motion delay and hold time, temperature-dependent operation and air function,
- **The incubator is equipped with:**
 - **Temperature regulator (max. 70°C),**
 - **Main switch-on module includes:**
 - Main switch to select the operating time mode,
 - Three unit standby lights (green, red, yellow),
 - Slider for fresh air regulation,
 - **Overheating regulator,**
- Additional glass door,
- **There are minor scuffs on the outer door of the incubator, visible in the photos,**
- The interior of the chamber is made of stainless steel,
- **Specification:**
 - Chamber's capacity: 749 litres,
 - Maximum shelves number: 10
 - Maximum one shelf load: 30 kg,
 - Maximum overall load: 160 kg,
- A mobile incubator on 4 wheels, two of which have breaks,
- **The set includes:**
 - Net trays; dimensions: 56 x 104 – 3 pieces,
 - English instruction manual (PDF),
- Dimensions:
 - The chamber: 60 x 103,5 x 120 cm,
 - The incubator: 76 x 119 x 162,5 cm,
- Weight: 187,5 kg,
- Provided with a valid Technical Passport,
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,



