



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STIEGELMEYER 6074-03-85-00-0 bedside cabinet (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Germany,
- **Stiegmeyer** is a stylish bedside cabinet that can be used in medical centers, hospitals, etc.,
- Ideal for rooms with limited space,
- Made of particle board with veneer finish,
- Wooden components made of laminated board,
- Double-sided opening drawers and door allowing the cabinet to be placed on both sides of the bed,
- Easy to keep it clean,
- Rounded edges on all sides of the cabinet, protect against damage,
- **It has:**
 - The table top with rounded edges with the following dimensions: 49.5 x 49.5 cm,
 - A drawer equipped with a pull-out separator with the following dimensions: 40 x 37.5 x 8.5 cm,
 - An open shelf with the following dimensions: 44 x 40 x 20,5 cm,
 - A closing space divided into two levels - with a sliding shelf inside, allowing the patient to adjust,
 - A basket on the door,
- A mobile cabine, on 4 plastic wheels,
- Overall dimensions: 50 x 50 x 90 cm,
- Weight: 31 kg,
- Cat 02,
- Provided with a valid Technical Passport,
- Guarantee:
 - 6 months for the domestic market (Poland),
 - 3 months for the international market,
- Financing possibilities (Poland only): Installments, leasing, loan,

