



### **KARL STORZ Model AUTOCON 350 surgical diathermy (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 300 W,
- **KARL STORZ AUTOCON 350** surgical diathermy is used in all areas of surgery. For safer effects on the patient's body, it provides soft and standard coagulation, as well as forced coagulation. It is characterized by precise power control. The device operates in monopolar cutting, monopolar and bipolar coagulation modes,
- **In the set:**
  - KARL STORZ foot attachment,
  - Brand new monopolar electrode cable with BOWA ERGO PEN holder - 2 buttons, 2.4 mm shaft, 4.5 m long (no electrodes included, can be extra delivered at additional cost),
  - Brand new BOWA clip cable for disposable neutral electrodes, 4.5 m long (no electrodes included, can be extra delivered at additional cost),
  - Brand new bipolar cable, BOWA tweezers 4 m long (bipolar tweezers can be delivered at additional cost),
  - Power cord,
- Dimensions: 45 x 33 x 17 cm,
- Weight: 11 kg,
- Provided with a valid Technical Passport,
- Guarantee:
  - 6 months for the domestic market (Poland),
  - 3 months for the international market,
- Financing possibilities (Poland only): Installments, leasing, loan,

