



JOUAN EB 170 laboratory incubator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in France,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 300 W,
- **EB 170** is a laboratory incubator with forced convection,
- **Specification:**
 - Smooth temperature regulation ranging between: 30o - 60o C,
 - Smooth ventilation regulation in the range of: 1-10,
 - The chamber made of stainless steel,
 - The chamber's capacity: 170 L,
 - Extra equipped with glass door,
 - Maximum number of shelves: 13,
- **In the set:**
 - Sieve trays with the following dimensions: 47 x 47 - 2 szt.,
 - Power cord,
- Dimensions:
 - Of the incubator: 65 x 63,5 x 89 cm,
 - Of the chamber: 51 x 51 x 62 cm,
- Weight: 70 kg,
- Provided with a valid Technical Passport,
- Guarantee:
 - 6 months for the domestic market (Poland),
 - 3 months for the international market,
- Financing possibilities (Poland only): Installments, leasing, loan,

