



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



DRAGER MIB II PROTO CONV-LS 7256931 protocol converter (digital adapter) (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Germany,
- Power: 2 W,
- The converter allows Drager's Infinity Series patient monitors to be joined with other devices in healthcare units - such as ventilators, incubators and anesthesia systems. Data from these devices is then displayed on the monitor, where it can be viewed and used to create trends and calculations,
- Economic solution for data integration,
- **Dedicated to use together with the following patients' Infinity monitors:**
 - Delta,
 - Delta XL,
 - Kappa,
- **It has in-built:**
 - RS 232 input (25 - pin subminiature D-type connector),
 - RJ45 MIB output,
- Dimensions: 8,4 x 15 x 2,8 cm,
- Weight: 0,4 kg,
- Provided with a valid Technical Passport,
- Guarantee:
 - 6 months for the domestic market (Poland),
 - 3 months for the international market,
- Financing possibilities (Poland only): Installments, leasing, loan,

