



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STRYKER Model VISION ELECT HDTV 26" medical monitor (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the USA ,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 150 W,
- **STRYKER Vision Elect HDTV** is an intelligent, microprocessor-based TFT-LCD monitor designed for the use in endoscopic surgical applications. It has an ergonomically designed display and is compatible with most RGB display standards,
- **Specification:**
 - LCD display,
 - Diameter: 26",
 - Color's depth: 16 millions,
 - Aspect ratio: 16:9,
 - Brightness: 500 cd/m²,
 - Reaction time: <25ms,
 - Viewing angle: +/- 89°,
 - Resolution: 1920 x 1200 @ 60Hz,
 - Input signal: DVI, VGA, HD/SD-SDI, C-Video/SOG, S-Video, Component (Y/G, Pb/B, Pr/R, H/CS, VS),
 - English menu,
- VESA type mounting (handle and stand are not included),
- **In the set:**
 - English instruction manual (PDF),
 - 24 V power adapter,
 - Power cord,
- Size or overall dimensions: 62 x 43 x 12 cm,
- Weight: 13 kg,
- Provided with a valid Technical Passport,
- Guarantee:
 - 6 months for the domestic market (Poland),
 - 3 months for the international market,
- Financing possibilities (Poland only): Installments, leasing, loan,

