



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **FAMED ŻYWIEC WS-29.5 urological attachment (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made Poland,
- **It comprises of:**
  - A bowl made of stainless steel, with the following dimensions: 65 mm x 265 mm x 325 mm and capacity of 4 liters,
  - Metal base made entirely of stainless steel with a matte finish. The base of the attachment allows it to be mounted to the operating table and rotated around the vertical axis,
- Overall dimensions: 27 x 41 x 26 cm,
- Set's weight: 4 kg,
- Provided with a 12-month valid Technical Passport,
- Guarantee:
  - 6 months for the domestic market (Poland),
  - 3 months for the international market,
- Financing possibilities (Poland only): Installments, leasing, loan,

