



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



THERMO SCIENTIFIC JOUAN INNOVENS 234 EB1 CAB laboratory incubator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in France,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 270 W,
- **EB1 Innovens** is a laboratory incubator with natural convection,
- Triple isolation and carefully placed heating elements ensure excellent temperature uniformity throughout the chamber,
- **Specification:**
 - Smooth temperature regulation in the range of: 5o - 75o C,
 - The chamber made of stainless steel,
 - Chamber's capacity: 234 L,
 - Equipped with extra glass door,
 - Maximum number of shelves: 14,
 - Microprocessor controls,
 - Digital temperature display,
- **In the set:**
 - Sieve trays with the following dimensions: 45 x 45 - 4 szt.,
 - Power cord,
- **Dimensions:**
 - Of the incubator: 68 x 63,5 x 120 cm,
 - Of the chamber: 50 x 50 x 94 cm,
- Weight: 80 kg,
- Provided with a valid Technical Passport
- **Guarantee:**
 - 6 months for the domestic market (Poland),
 - 3 months for the international market,
- Possible guarantee extension for extra fee, up to 12 months for an international market,
- Financing possibilities (Poland only): Installments, leasing, loan,

