



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



STOLTER U-03 handgrip (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Poland,
- **Stolter U-03** is a hand support handle designed for universal general examination tables and gynecological chairs,
- **In the set:**
 - Mounting supports for the rail from the table – 2 pieces,
- Overall dimensions: 8,5 x 19,5 x 36 cm,
- The set's weight: 3 kg,
- Provided with a valid Technical Passport
- Guarantee:
 - 6 months for the domestic market (Poland),
 - 3 months for the international market,
- Possible guarantee extension for extra fee, up to 12 months for an international market,
- Financing possibilities (Poland only): Installments, leasing, loan,

