



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ZOLL M SERIES two-phase defibrillator (Reconditioned)

- Recodnitioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the USA,
- Power supply: 10 – 29 V and via 1, new in-built 10V / 2,5 Ah battery,
- Frequency: 50/60 Hz,
- Power: 130 W,
- **Zoll M Series** are lightweight and portable devices that combine the functions of defibrillation, ECG, patient's monitoring with the ability to print and record data,
- The defibrillator is designed for use during any type of resuscitation, and its light weight, small size and compact design make it particularly suitable for the use during transport,
- **It has:**
 - 5,5" LCD display,
 - English menu,
 - English keyboard on the panel control,
 - **PACER (STIMULATOR) function - switches on the stimulation function,**
 - **The following modules and wirings are in-built:**
 - SPO2 – clip sensor,
 - ECG – 10-lead wire ECG,
 - RS-232/EKG X 1000 connector,
- A printer for paper of 90 mm width,
- PCMCIA external memory cards port,
- Portable thanks to convenient handle and a bag for wires,
- **In the set:**
 - Dual-purpose (adult-child) hard paddles,
 - Polish instruction manual (PDF),
 - Brand new power cord,
- Overall dimensions: 36 x 47 x 22 cm,
- The set's weight: 10 kg,
- Cat D121,
- Provided with a valid Technical Passport,
- Guarantee:
 - 6 months for the domestic market (Poland),
 - 3 months for the international market,
- Financing possibilities (Poland only): Installments, leasing, loan,

