



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ZOLL M SERIES two-phase defibrillator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the USA,
- Power supply: 10 – 29 V and via 1 in-built 10 V / 2,5 Ah battery,
- Frequency: 50/60 Hz,
- Power: 130 W,
- **Zoll M Series** are lightweight and portable devices that combine the functions of defibrillation, ECG, patient's monitoring with the ability to print and record data,
- The defibrillator is designed for use during any type of resuscitation, and its light weight, small size and compact design make it particularly suitable for use during transport,
- **It has:**
 - 5,5" LCD display,
 - English menu,
 - English keyboard on the control panel,
 - **An in-built module:**
 - ECG – 4-lead ECG wire,
 - RS-232/EKG X 1000 connector,
 - A printer for paper of 90 mm width,
 - External memory cards port PCMCIA,
- Portable thanks to a convenient handle and a bag for wires (damaged pocket zipper - visible in the photo),
- **In the set:**
 - Dual-purpose (adult-child) hard paddles,
 - Polish instruction manual (PDF),
 - Brand new power cord,
- Overall dimensions: 36 x 47 x 22 cm,
- Set's weight: 10 kg,
- Cat D122,
- Has a valid Technical Passport,
- Guarantee:
 - 6 months for the domestic market - Poland,
 - Months for the international market,
- Financing possibilities (Poland only): Installments, leasing, loan,

