



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



O2 MEDAP FINA FLOW C2 15 single oxygen dispenser (Ex-Display)

- A new, post-display product,
- **Real photos of the product,**
- Made in Germany,
- **Fina Flow oxygen dispenser** with pressure compensation is designed to regulate the flow of properly humidified medical oxygen supplied from the central gas system to the patient,
- Connected via a dedicated plug-in tip to suitable medical gas outlet points (mounted directly in the wall, bedside panels, anesthesia columns, etc.),
- Large control valve allows extremely precise dosing during inhalation or blowing,
- Corrosion-resistant, matte-chrome plated brass and stainless steel ensure long service life,
- **Specification:**
 - Gas type: oxygen,
 - Flow range: 0 - 15 l/min,
 - Accuracy:
 - $\pm 10\%$ above 5 l/min flow rate indicator,
 - ± 0.5 l/min up to 5 l/min flow rate indicator,
 - Connection system: wall-mounted MEDAP,
 - Connection threads: UNF 9/16",
 - Operating pressure: 4,5 - 5,5 bar,
- **In the set:**
 - Brand new bubble humidifier with safety valve,
- Dimensions:
 - Of the dispenser: 14 x 5 x 17 cm,
 - Overall: 17 x 7,5 x 34 cm,
- Set's weight: 830 g,
- Provided with a valid Technical Passport,
- Guarantee:
 - 6 months for the domestic market (Poland),
 - 3 months for the international market,
- Financing possibilities (Poland only): Installments, leasing, loan,

