



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



FROILABO FIRLABO BIO BCS245 incubator (Reconditioned)

- Recodnitioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in France,
- Power supply: 230 V,
- Frequency: 50-60 Hz,
- Power: 500 W,
- **FROILABO FIRLABO BIO BCS245** is a biological incubator with natural Convection,
- **Most important advantages:**
 - Temperature range: from +5°C to 65°C of ambient temperature,
 - High-performance PID electronic control,
 - Overheating protection according to the EN61010 standard,
 - Control by Pt 100 Ohm probe,
 - Digital display with 0.1°C accuracy (real value and setpoint),
 - Stainless steel chamber interior with rounded edges for easy cleaning and disinfection,
 - Glassed Plexiglas interior door to reduce temperature drops,
 - Height-regulated shelves,
- **Specification:**
 - Capacity: 245 l,
 - Nominal temperature: od +5°C do 65°C,
 - Environmental uniformity at 37°C (±°C): <1,
 - Environmental uniformity at 37°C (±°C): <1,65°C (±°C): <2,
 - Time stability (±°C): <0,2,
 - Shelves regulation levels: 62,
 - Maximum load per shelf: 20 kg,
- **In the set:**
 - Stainless steel perforated shelves (3 pcs.),
 - Power cord,
- Dimensions:
 - Overall: 64 x 62 x 124 cm,
 - Of the chamber: 50 x 50 x 98 cm,
 - Of the shelf: 45 x 46 cm,
- Maximum load: 90 kg,
- Weight: 80 kg,
- Provided with valid Technical Passport,
- Guarantee:
 - 6 months for the domestic market (Poland),
 - 3 months for the international market,
- Financing possibilities (Poland only): Installments, leasing, loan,

