



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDTRONIC BIPHASIC LIFEPAK 12 defibrillator (Reconditioned)

- Recodnitioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the USA,
- Power supply: 230 V or via 2, exchangeable 12V 1,2 Ah batteries (**one included in the set, possibility to extra equip with the second one at additional fee**),
- Frequency: 50-60 Hz,
- Power: 180 W,
- **Lifepak 12** is a two-phase ambulance defibrillator,
- It has an 8.1" display,
- **Menu to be chosen in the following languages: PL, EN, FR, DE, ES, SV, IT, FI, NL, PT, DA, NO,**
- English keyboard on the control panel,
- An in-built printer for 11 cm paper width,
- Equipped with a transportation handle,
- **Functions:**
 - SYNC - Synchronization - Cardioversion,
 - ADVISORY - INFO Enables the Continuous Patient Surveillance System (CPSS),
 - ANALYZE - Activates the SAS (Shock Advisory System),
 - PACER - Stimulator - Activates the stimulation function,
- **Equipped with the following modules with wirings and sensors:**
 - ECG - 5 - lead cable,
 - SPO₂ - brand new sensor clip,
- **In the set:**
 - Hard paddles,
 - Polish instruction manual (PDF),
 - Medtronic charger - power supply 230 V,
 - Power cord for a charger,
- Overall dimensions: 23 x 40 x 38 cm,
- Set;s weight: 10 kg,
- Cat. 65,
- Provided with valid Technical Passport,
- Guarantee:
 - 6 months for the domestic market (Poland),
 - 3 months for the international market,
- Financing possibilities (Poland only): Installments, leasing, loan,

