



NIDEK KM-300 keratometer ophthalmometer (Reconditioned)

- Recodnitioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Japan,
- Power supply: 230 V,
- Frequency: 50-60 Hz,
- Power: 20 W,
- The **NIDEK KM-300 keratometer** is used to measure the curvature of the frontal surface of the cornea which helps in assessing the degree and axis of astigmatism,
- Fully manual,
- Massive and stable construction ensures precision measurement,
- Dimensions: 48 x 18 x 40 cm,
- Weight: 9,5 kg,
- Provided with valid Technical Passport,
- Guarantee:
 - 6 months for the domestic market (Poland),
 - 3 months for the international market,
- Financing possibilities (Poland only): Installments, leasing, loan,

