



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MEDTRONIC BIPHASIC LIFEPAK 12 defibrillator (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the USA,
- Power supply: 12 V and via two exchangeable 12V 1,2 Ah batteries (**one is included, the second one can be added for extra fee**),
- Frequency: 50-60 Hz,
- Power: 180 W,
- **Lifepak 12** is a two-phase ambulance defibrillator,
- It has an 8.1" display,
- **Menu to be chosen in the following languages: PL, EN, FR, DE, ES, SV, IT, FI, NL, PT, DA, NO,**
- English keyboard on the control panel,
- An in-built printer for 11 cm paper width,
- Equipped with a transportation handle,
- **Functions:**
  - SYNC - Synchronization - Cardioversion,
  - ADVISORY - INFO - Enables the Continuous Patient Surveillance System (CPSS),
  - ANALYZE - Activates the SAS (Shock Advisory System),
  - PACER - Stimulator - Activates the stimulation function,
- **Equipped with the following modules with wirings and sensors:**
  - ECG - 10-lead cable,
- **In the set:**
  - Hard paddles,
  - Polish instruction manual (PDF),
  - Medtronic charger – power supply of 230 V,
  - Power cord for the charger,
- Overall dimensions: 23 x 40 x 38 cm,
- Set's weight: 10 kg,
- Cat. 66,
- **It has up-to-date overview and is ready to work,**
- **Provided with a 12-month valid Technical Passport (service report),**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months for the international market,**
  - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

