



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS ENF-P4 rhinolaryngoscope (Reconditioned)

- Reconditioned (used),
- Technical condition: good,
- Visual condition: good,
- **Real photos of the product,**
- Made in Japan,
- **ENF-P4** - provides clear, comprehensive images of the nasal cavity and laryngeal and pharyngeal areas,
- Uses high-quality fiber optic technology to provide clear and accurate images, while allowing us to add even more fibers for larger and brighter images,
- The device has blind spots **visible in the photo,**
- **Specification:**
 - Viewing direction: forward view,
 - Maximum angle of field of view: 85°,
 - Outer diameter of tip: 3.6 mm,
 - Working length: 300 mm,
 - Overall length: 530 mm,
 - Maximum endoscope tip bending angle:
 - Up: 130°,
 - Down: 130°,
- **In the set:**
 - Olympus case,
- Case's dimensions: 9 x 42 x 32 cm,
- Set's weight: 2 kg,
- Cat 02,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport and ready to work,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,



