



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



THREMO SCIENTIFIC FORMA FFGL430V19 laboratory freezer (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the USA,
- Power supply: 230 V,
- Emergency supply: via 1 new in-built 12V 1.2 Ah battery,
- Frequency: 50/60 Hz,
- **THREMO SCIENTIFIC FORMA** is a high-performance laboratory freezer with automatic defrost function, designed to operate at -20 °C,
- **Specification:**
 - Powder-coated housing,
 - CFC-free, high-density polyurethane foam insulation protects the enclosure from ambient temperature fluctuations and minimizes operating costs,
 - Durable door with key lock for extra security (no key included),
 - Self-closing door with magnetic seal,
 - Epoxy resin coated wire shelves regulated in 2.5 cm increments,
 - Automatic defrost,
 - Digital temperature display,
 - An in-built alarm (audible and visual warnings),
 - Access port: 25 mm,
 - Capacity: 133 l,
 - Has one operating program with min. operating temperature: - 20 °C,
 - Mobile on 4 wheels with a brake,
- **English keyboard on the control panel,**
- **In the set:**
 - Power cord,
 - 2 x wire shelves
 - Key to the switch,
- Dimensions:
 - Overall: 68 x 61 x 90 cm,
 - Of the chamber: 50 x 50 x 40 cm,
 - Of the shelf: 48 x 49 cm,
- Weight: 80,5 kg,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport and ready to work,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

