



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



VERMEIREN EAGLE 620 patient lift for bathing and transportation (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Belgium (2018),
- Battery power supply: via 2 new, in-built 12V 2,9 Ah batteries,
- Frequency: 50/60 Hz,
- **EAGLE 620 VERMEIREN patient lift** is used for easy transport and lifting of patients. The lift proves useful when transferring a patient, for example: from a bed to a wheelchair. It also helps with hygiene and nursing procedures, bathing and other activities related to moving and transporting patients. The undoubted advantages of the EAGLE lift are quick operation and easy folding for transport or storage. An additional advantage is a high-quality, durable actuator that guarantees smooth and trouble-free operation,
- **Specification:**
 - Durable design that allows lifting capacity up to 150 kg,
 - Includes REGULAR size L (70-120 kg) carrier (seat),
 - The lift works automatically controlled by remote control,
 - High-quality LINAK actuators guarantee smooth and trouble-free operation,
 - The design allows easy and quick folding for transport without using any tools,
 - The sliding mechanism of a frame is movable and allows smooth regulation of spacing for maximum stability,
 - Lift equipped with emergency button,
 - The design of the lifts allows patients to be lifted, e.g. from a bed, from a wheelchair,
 - Electric arm height regulation in the range of: 72 - 174 cm,
 - Distance between legs manually regulated in the range of: 55 - 103 cm,
 - Number of operations: about 40 lifts per load,
 - Equipped with 4 wheels, 2 of which have a locking mechanism,
- **In the set:**
 - Instruction manual in the following languages: PL, EN, FR, NL, DE, IT, ES (PDF),
 - Charger,
 - REGULAR Carrier, Size L (70-120kg)
- Dimensions:
 - Base length: 111 cm,
 - Base height: 110 cm,
- Maximum load: 150 kg,
- Weight: 38 kg
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport and ready to work,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

