



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



BQE DYNAMIC doctor's chair (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the Netherlands,
- **Dynamic** is an ergonomic medical chair that activates and corrects the sitting position. It features a unique patented scissor mechanism between the seat and backrest. The lumbar support automatically moves backward and adjusts to the empty contours. The natural S-shaped anatomical shape of the spine remains even when sitting passively. The small lumbar support provides optimal freedom of movement and activates spinal stretching,
- **Durable, medically certified upholstery in black, has small traces of use (there is a possibility to change upholstery to new one from our color pattern book, for extra fee),**
- Manual seat height regulation by gas spring in the range of: 55-71 cm,
- Perfect for doctors' offices and laboratories or, for example, as a dental chair,
- The base of the cross is made of chrome-plated steel, equipped with five double wheels,
- **In the set:**
 - Instruction manual in the following languages: EN, DE, NL,
- Dimensions:
 - Of the seat: 37,5 x 43,5 x 9 cm,
 - Overall: 60 x 60 x 76 cm,
- Maximum load: 110 kg,
- Weight: 11,5 kg,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (service report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

