



ALVO 4-09-100 operating table leg holder (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Poland,
- **ALVO orthopedic attachment** is used to operate on patients with fractures of lower limb bones. The use of additional equipment gives the possibility to perform specialized procedures. Metal elements are made of stainless, acid-resistant steel. Changing the length and angle of opening of the arms of the attachment is possible thanks to manually adjustable locking mechanisms. The arrangement of the arms of the attachment provides the possibility of monitoring the patient's lower limbs. The appropriately shaped top of the attachment is mounted in place of the footrests in the top of the operating table. It has the shape of a triangle with a cut apex. The location of the stretching appliance and the resistance roller ensures effective stretching of the operated limb. The tension appliance itself, thanks to its functions, ensures precise positioning and stretching of the leg during surgery. The design of the attachment allows limb monitoring with the C-arm.
- Easy transportation and space saving due to the ability to store all accessories on the cart,
- **Medical-approved upholstery in black (upholstery can be replaced with new one, at additional cost),**
- Dimensions: 215 x 60 x 110 cm,
- Weight: 63 kg,
- **It has an up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

