



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



J&J CIDEX 2027 sterilization container (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the USA,
- **J&J Cidex** is a container for sterilization and disinfection of medical instruments,
- Can be used for manual cleaning or high level disinfection,
- Temperature resistant up to 121°C,
- Reduces the evaporation of steam and vapors due to the unique rounded edge,
- Easy cleaning thanks to the wide radius corner
- Capacity: 14.1 L,
- Blue,
- **In the set:**
 - Perforated inner tray with dimensions: 15 x 63 x 7 cm,
 - The cover,
- Dimensions:
 - Internal: 14 x 66,5 x 13,4 cm,
 - Overall: 19 x 72,3 x 13,8
- Weight: 2,5 kg,
- Cat 02,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport and ready to work,**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

