



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PHYSIO-CONTROL (3010267-00) top pouch for LIFEPAK 12/15 defibrillator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the USA,
- A pouch designed to store sensors and electrodes for the Lifepak 12 and 15 defibrillator,
- Attached instead of standard hard paddles,
- Equipped with a zipper,
- Dimensions: 34 x 11 x 8,5 cm,
- Weight: 0,2 kg,
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

