



ELKON CSL-75 laboratory incubator - dryer (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Poland,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 0,800 kW,
- **An in-built microprocessor regulator LUMEL RE3 controller:**
 - Process value display,
 - Setpoint display,
 - Alarm light,
- Maximum temperature: 110° C,
- **The chamber equipped with:**
 - Protective glass,
 - 4 pcs. - tray holders, with the ability to hang them in 12 height levels,
 - 2 pcs. - mesh trays with dimensions: 35.5 x 43 cm,
- Dimensions:
 - Of the chamber: 37 x 43 x 37 cm,
 - Overall: 62 x 71 x 56,5 cm,
- Weight: 64,5 kg,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (service report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

