



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **MEDTRONIC BATTERY SUPPORT SYSTEM 2 docking station - charger (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the USA,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 120 V,
- **MEDTRONIC Battery Support System 2** is a highly intelligent and automated system that recognizes and optimizes the maintenance of NiCd batteries (FASTPAK, FASTPAK 2 and Lifepak NiCd), as well as sealed lead acid batteries (Lifepak SLA). It includes recognition software to charge SLA batteries along with NiCd batteries. All three battery chambers perform charging, conditioning and battery life tests.
- **Specification:**
  - Optimal battery performance provided by improved NiCd pulse charging algorithm,
  - NiCd battery conditioning cycles ensure long battery life,
  - User-friendly system notifies users when batteries are ready for use or need to be discarded,
  - Ability to add new features and enhancements as future battery technology evolves,
- **English keyboard on the control panel,**
- **In the set:**
  - English instruction manual (PDF),
  - Power cord,
- Dimensions: 51 x 41 x 31 cm,
- Weight: 8 kg,
- Cat 02,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (service report),**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months for the international market,**
  - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

