



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



One-sided X-RAY protection apron with a STAR X-RAY collar (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the USA,
- **Star X-RAY protection apron** for radiation protection is designed to combine comfort and uncompromising secondary radiation protection,
- Made of lightweight multilayer and flexible lead rubber for superior performance.
- Protects the front, shoulders and thyroid of the user, mainly during dental x-rays,
- Fastens with Velcro,
- Radiation attenuation equivalent: front 0.30 mm Pb,
- Dimensions:
 - Width: 61 cm
 - Length: 90 cm,
 - Overall (when folded) 31 x 23,5 x 4 cm,
- Weight: 2,5 kg,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (service report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

