



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



A-DEC 1005 dental chair (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the USA,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- **The A-DEC 1005 is an electrically operated hydraulic dental chair with very quiet operation,**
- **Medically certified upholstery in gray (it will be replaced with new upholstery in the color of the customer's choice from our color pattern book),**
- All structural components are metal, which ensures durability and firmness,
- The base is made of thick solid steel plate for strength and stability,
- **Electrically controlled with a foot switch and built-in buttons in the seat back:**
 - Seat height in the range of: 45 - 82 cm,
 - Backrest angle,
- **Manual regulation of:**
 - Headrest,
 - Armrests that slide backwards to make it easier for patients to enter and exit from either side,
 - Swivel chair base,
- Dimensions:
 - Of the seat: 48 x 120 cm,
 - Of the backrest: 48 x 50 cm,
 - Of the base: 71 x 84 cm,
 - Overall: 76 x 130 x 112 cm,
- Weight: 142 kg.
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (service report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

