



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MELAG VACUKLAV 31 type autoclave (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 2500 W,
- **MELAG Vacuklav 31** is an efficient steam autoclave offering 4 operating programs and 2 test programs,
- **Product's characteristics:**
 - B class,
 - Automatic,
 - Pre-vacuum and fractional vacuum generated by a vacuum pump,
 - Additional vacuum drying is possible,
 - Preheat option helps to maintain temperature between cycles,
 - Equipped with indicator lights to indicate unit errors,
 - 2.5" LCD display,
 - Foiled control panel,
 - On the lid of the autoclave, a place to put down the tray,
 - Water drain connector on the front of the unit,
 - Sockets:
 - For Melaprint printer (not included),
 - For Meladem ion exchange system (not included),
 - 230 V outlet,
 - Made of stainless steel materials,
 - Stable on support legs, including front legs with height regulation,
 - German keyboard on the control panel,
- **It has 4 operating programs and 2 test programs:**
 - Unpacked tools: pressure 2 bar - time 4 min - temperature 134 °C,
 - Packed tools: pressure 2 bar - time 5 min - temperature 134 °C,
 - Textiles: pressure 1 bar - time 17 min - temperature 121 °C,
 - Disinfection: pressure 0.3 bar - time 8 min - temperature 105 °C,
 - Bowie & Dick test: vapor permeation test to check the efficiency of the autoclave by the user,
 - Leakage test: vacuum test to check the leakage of the autoclave,
- **Specification:**
 - Sterilization chamber capacity: 14,5 L,
 - Maximum load capacity: 4 kg of instruments including trays or 1.5 kg of textiles,
 - Water consumption during sterilization: 300 - 400 ml, depending on the program and loading level,
- **In the set:**
 - English instruction manual (PDF),
 - Trays of dimensions: 18 x 28 x 2 cm - 2 pieces,
 - Trays' holder,
 - A stand with holders for 3 trays,
 - Brand new condensate tank,
- Dimensions:
 - Of the autoclave: 63 x 45 x 45 cm,
 - Of the chamber: Ø 23 cm x 34 cm,
- Weight: 44 kg,
- Cat 04,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (service report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

