



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



SONICAID FM830 ADAPT CTG device (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in England,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 100 VA,
- **SONICAID FM800 series cardiotocograph** is a high-quality cardiotocograph suitable for monitoring single or twin pregnancies with simultaneous monitoring of maternal parameters. In addition, the device is equipped with two analyses of CTG recording (antenatal Care and postnatal Trnd),
- **Specification:**
 - In-built 320 x 240 TFT Liquid Crystal Display,
 - Internal memory,
 - In-built alarms,
 - Function of automatic movement detection and determination of fetal movement activity,
 - Can be connected to digital central monitoring systems, telemetry and external computer software,
 - In-built thermal printer,
 - Ability to monitor twins simultaneously,
- **Menu to be chosen in: EN, ES, IT, DK, FR, DE, NL**
- **English keyboard on the control panel,**
- **In the set:**
 - English instruction manual (PDF),
 - Brand new toco head,
 - Brand new US head,
 - NBIP cable + cuff (25-35 cm),
 - Brand new temperature sensor,
 - Systolic patient button,
 - Power cord,
- Dimensions: 37 x 39 x 20 cm,
- Weight: 10 kg,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

