



## **WELCH ALLYN PROBP 2400 digital blood pressure device (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the USA,
- Power supply: 230 V,
- Emergency power supply: via 4 new, in-built AA 1,2V 2400 mAh batteries,
- Frequency: 50/60 Hz,
- **WELCH ALLYN PROBP 2400** is a digital blood pressure monitor that provides fast, accurate and reliable digital readings at an affordable price. Ideal for primary care and specialized outpatient settings, the PROBP 2400 is equipped with blood pressure averaging technology that provides a more complete picture when detecting and diagnosing hypertension,
- **Specification:**
  - Detection of irregular heartbeat provides data to help with decision-making for testing,
  - Microprocessor-controlled deflation rate provides consistent technique,
  - Averaging mode of three readings can give a broader picture of a patient's actual blood pressure,
  - Manual mode allows the ProBP 2400 to serve as a digital meter when auscultation is indicated,
  - Lithium-ion battery technology provides up to 100 readings on a single charge,
  - BP accuracy +/- 3mmHg. Pulse +/- 5%,
- Blood pressure device installed on a mobile stand with a 5-arm mobile base with wheels (all with brakes), a basket on the column for convenient storage of cuffs and tubing,
- **In the set:**
  - English instruction manual (PDF),
  - NBIP cable + cuffs: 25-34 cm, 20-26 cm, 15-21 cm,
  - Power supply,
- Dimensions: 51 x 51 x 106 cm,
- Weight: 6 kg,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months for the international market,**
  - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

