



BLACK OPTIC LTD RONDO ophthalmic table (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 210 VA,
- **Rondo** is an ophthalmic electric table designed to work with diagnostic ophthalmic equipment,
- **Equipped with:**
 - Durable straight work surface with dimensions: 50 x 90 cm (**has minor signs of use**),
 - Two buttons on the front of the table:
 - On/Off switch,
 - To regulate the height of the tabletop in the range: 70 - 110 cm,
 - Two 230V outlets, located under the tabletop, to power table-top appliances,
 - Double wheels with brakes,
- Overall dimensions: 46 x 66 x 70/89 cm,
- Weight: 37,5 kg,
- Cat 01,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

