



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



LEITZ WETZLAR LABORLUX K laboratory microscope (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Germany,
- Power supply: 230 V,
- **LEITZ WETZLAR LABORLUX K** laboratory microscope,
- Smooth regulation of the intensity and focus of the light,
- Light source: 6V 20W halogen bulb,
- Two LEITZ PERIPLAN 10x/18 eyepieces with smooth regulation of spacing,
- A 5-slot revolver with three PHACO lenses:
 - 10/0.25,
 - 40/0.65,
 - 100/1.25-0.60,
- Scale table with paddles for smooth movement in X/Y surfaces using two knobs,
- Smooth height regulation of the condenser,
- **In the set:**
 - Power cord,
- Dimensions: 33 x 21 x 39 cm,
- Weight: 9 kg,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

