



IONTO-STERIL UL dry air sterilizer (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Germany (2015),
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 690 W,
- **IONTO-STERIL UL** is a hot air sterilizer, it is a fully safe and durable device. It meets all the requirements for compact devices,
- Safety is guaranteed by the automatic shutdown function of the sterilization cycle including power off when the cycle is completed,
- The device can be used in veterinary and medical applications,
- Dry sterilization,
- Forced air circulation,
- Capacity: 4.3 L,
- Smooth temperature control of up to max. 180° C,
- Has a timer switch adjustable from: 0 - 120 min,
- Stainless steel case, strong construction and long life of the device,
- **In the set:**
 - 2 x trays,
 - Tray holder,
 - Power cord,
- Dimensions:
 - Overall: 31 x 41,5 x 26 cm,
 - Of the tray: 18 x 28 cm x 2 cm,
- Weight: 14 kg,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

