



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **JETPEEL-3 a high pressure oxygen delivery device (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Israel,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 280 W,
- Using a unique cap with a built-in nozzle, **JETPEEL-3** produces a jet of gas and liquid that interacts with the skin. At the point of contact with the skin, a deformation in the shape of a cavity is created and erosion is produced at the bottom of the cavity. In this way, a peeling effect is obtained. Static pressure and coolness around the jet provide a pain-relieving effect. While the procedure is being performed, the patient only feels a pleasant wrapping with a strong stream of cool air.
- **Oxygen micro-dermabrasion:**
  - Allows peeling of any depth on any skin,
  - No side effects or acute inflammation,
- **The system is used:**
  - For skin exfoliation to reduce surface and deep wrinkles,
  - For correcting scars in combination with other cosmetic treatments,
  - For acne ("acne") and acne lesions ("post-acne"),
  - For superficial and deep peeling,
  - For massage for lymphatic drainage and soft tissue lifting of the face and neck,
  - For scalp treatments,
  - For local skin moisturization,
  - For fast and effective skin cleansing,
  - For effective smoothing of scars and fighting cellulite,
  - For fighting stretch marks,
  - In post-surgical rehabilitation,
- **English menu,**
- **In the set:**
  - Foot attachment,
  - New treatment nozzle + hose,
  - New O2 hose,
  - Power cord,
  - Transport box,
  - **No reducer - can be added at extra cost,**
- Dimensions: 20 x 40 x 68 cm,
- Weight: 35 kg,
- Cat 01,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months for the international market,**
  - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

