



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OLYMPUS CYF-V2 video cystofiberscope (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Japan,
- **Specification:**
 - Viewing direction: forward view,
 - Outer diameter of the tip: 5.4 mm,
 - Operating channel diameter: 2.2 mm,
 - Operating length: 380 mm,
 - Overall length: 650 mm,
 - Depth of acuity: 3-50 mm,
 - Max. bend angle of the endoscope tip up and down: 210°/120°
 - Max. angle of field of view: 120°,
- **In the set:**
 - Case,
 - Valves,
 - English instruction manual,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

