



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **OLYMPUS WM-NP1 endoscopic cart (Reconditioned)**

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the UK,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- **Olympus WM-NP1** is a stable and durable cart suitable for carrying endoscopes and fragile medical equipment,
- **It has:**
  - Upper top with dimensions: 47 x 38 cm,
  - Lower top with dimensions: 47.5 x 56 cm,
  - 2 shelves, 40 x 54 cm, with manual height regulation (every 1.5 - 2 cm) with Allen screws,
  - Handles for easy movement of the cart,
  - Endoscope holder with manual height regulation,
  - Monitor holder with height regulation within 20 cm (6 points every 3.5 cm),
  - Electrical installation, 12 power sockets for the medical equipment transported,
- **Product's characteristics:**
  - Power switch on the front of the cart,
  - Detachable rear panel of the cart facilitates cable management,
  - Power cables aesthetically hidden in panels on the sides of the cart,
  - Corners of the shelves protected with safety pads, and edges over the wheels with tires,
  - Mobile on four wheels with Ø 12 cm, two of them with brakes,
- **In the set:**
  - Power cords,
- Overall dimensions: 67 x 66 x 166/200 cm,
- Weight: 86,5 kg,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months for the international market,**
  - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

