



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## OLYMPUS ENF-P4 rhinolaryngofiberscope (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Japan,
- **ENF-P4** - provides clear, comprehensive images of the nasal cavity as well as areas of the larynx and pharynx,
- Uses high-quality fiber-optic technology to provide clear and accurate images, while allowing us to add even more fibers for larger and brighter images,
- **The device has seven blind spots,**
- **Specification:**
  - Viewing direction: forward view,
  - Maximum field of view angle: 85°,
  - Tip outer diameter: 3.6 mm,
  - Operating length: 300 mm,
  - Overall length: 530 mm,
  - Maximum endoscope tip bending angle:
    - Up: 130°,
    - Down: 130°,
- **In the set:**
  - Olympus case,
- Case's dimensions: 9 x 42 x 32 cm,
- Set's weight: 2 kg,
- Cat 02,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months for the international market,**
  - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

