



OLYMPUS ENF-P4 rhinolaryngofiberscope (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Japan,
- **ENF-P4** - provides clear, comprehensive images of the nasal cavity as well as areas of the larynx and pharynx,
- Uses high-quality fiber-optic technology to provide clear and accurate images, while allowing us to add even more fibers for larger and brighter images,
- **The device has seven blind spots,**
- **Specification:**
 - Viewing direction: forward view,
 - Maximum field of view angle: 85°,
 - Tip outer diameter: 3.6 mm,
 - Operating length: 300 mm,
 - Overall length: 530 mm,
 - Maximum endoscope tip bending angle:
 - Up: 130°,
 - Down: 130°,
- **In the set:**
 - Olympus case,
- Case's dimensions: 9 x 42 x 32 cm,
- Set's weight: 2 kg,
- Cat 02,
- **It has up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

