



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## MELAG Type 17 autoclave (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 1310 W,
- Melag 17 is an N-class steam autoclave,
- **It has 2 operating programs:**
  - Operating pressure 1 bar - time 50 min - temperature 121° C,
  - Operating pressure 2 bar - time 30 min - temperature 134° C,
- On the lid of the autoclave, a place to put the tray away,
- Built-in pressure gauge and thermometer,
- Water drainage connector on the front of the unit,
- Pressure drain button,
- Indicator lights for autoclave start-up and heating,
- All steam guiding elements made of stainless materials,
- Height-regulated front feet,
- **Specification:**
  - Chamber capacity: 11 L,
  - Water consumption during sterilization: 300 ml,
  - Maximum load capacity: 3 kg of accessories including trays or 0.2 kg of textiles,
- English control panel,
- **In the set:**
  - Polish instruction manual (PDF),
  - Perforated tray, 40 x 14 x 2 cm (2 pcs),
  - Tray holder,
  - Drain hose,
  - Power cord,
- Dimensions:
  - Overall: 63 x 43,5 x 35,5 cm
  - Of the chamber: Ø 17,5 x 48 cm,
- Weight: 20 kg,
- Cat 02,
- **It has an up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months for the international market,**
  - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

