



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



ERBE IES 2 surgical smoke evacuator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Germany,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- **ERBE IES 2 surgical smoke evacuator** extracts harmful fugitive fractions from the surgical field and increases occupational safety by minimizing the amount of inhaled undesirable substances for medical personnel.
- The device cooperates with the VIO 200/300 series diathermy, enables effective smoke extraction from the operating field and its two-stage filtration. It is possible to direct the stream of filtered air to an external extraction system,
- Programmed and controlled by diathermy, but it is also possible to control the operation of the device manually through a control panel.
- Complete monopolar handpieces (disposable), Slim-line reusable handpieces and reusable hoses with replaceable tips are compatible with the device.
- **No main filter (can be added at extra cost),**
- **In the set:**
 - Power cord,
- Dimensions: 20,5 x 40 x 26,5 cm,
- Weight: 11 kg,
- **It has an up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

