



## DATEX-OHMEDA TEC 850 vaporizer (Isoflurane) (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the USA (2019),
- Designed to deliver isoflurane anesthetic, it combines significant performance, convenience and ergonomic benefits with important safety features that help improve delivery and reliability while reducing overall operating costs.
- **Main characteristics:**
  - For anesthesia applications - isoflurane application,
  - Precise administration according to the concentration controller over the entire clinical flow range.
  - Easy-to-use concentration controller and precise graduation ensure precise control of anesthetic application,
  - Optimized design for consistent and reliable performance over the life of the product,
  - TEC 850 vaporizer supports several popular plug-in adapter systems
  - Improved ergonomics and modern design as the ideal addition to GE Healthcare's anesthesia systems,
- **This particular vaporizer uses Selectatec (S - 2000) mounting type,**
- **Fill type: Easy Fill,**
- Dimensions: 11,5 x 22 x 26 cm,
- Weight: 7 kg,
- **It has an up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
  - **6 months for the domestic market (Poland),**
  - **3 months for the international market,**
  - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

