



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



MEDTRONIC BIPHASIC LIFEPAK 12 defibrillator (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in the USA,
- Power supply: 12 V and via two replaceable 12V 1,2 Ah batteries, **(one in the set, possibility to equip with the second one at extra fee),**
- Frequency: 50/60 Hz,
- Power: 180 W,
- **Lifepak 12** is a biphasic ambulance defibrillator,
- It has an 8.1" diagonal display,
- **Menu to be chosen in: PL, EN, FR, DE, ES, SV, IT, FI, NL, PT, DA, NO,**
- English keyboard on the control panel,
- Built-in printer for paper width: 11 cm,
- Equipped with a carrying handle,
- **Functions:**
 - SYNC - Synchronization - Cardioversion,
 - ADVISORY - INFO - Enables the Continuous Patient Surveillance System (CPSS),
 - ANALYZE - Activates the SAS (Shock Advisory System),
 - PACER - Stimulator - Activates the stimulation function,
- **Equipped with the following modules with wiring and sensors:**
 - ECG - 3-lead cable,
 - QUICK COMBO therapy cable,
- **In the set:**
 - Polish instruction manual (PDF),
- **The device can be extra equipped with the following (at additional cost):**
 - Hard paddles,
 - Extra battery,
 - External charger from AD Elektronik shown in the last photo,
- Overall dimensions: 22 x 39 x 32 cm,
- Set's weight: 7,5 kg,
- Cat. 70,
- **It has an up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

