



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



OPTOMIC OP-C2 (vision track) colposcope (Reconditioned)

- Reconditioned (used),
- Technical condition: very good,
- Visual condition: very good,
- **Real photos of the product,**
- Made in Spain,
- Power supply: 230 V,
- Frequency: 50/60 Hz,
- Power: 170 W,
- Light source: 15V 150W halogen bulb,
- **OP-C2 OPTOMIC optical colposcope** is designed for diagnostic observation, daily examinations in the gynecological office as well as small surgical centers. It features a very good ergonomic design, free head regulation and high-end Japanese optics. Thanks to its cooperation with a digital video track, you can observe a very high quality image on a TV monitor or computer. It is equipped with a series of precise mechanical systems for accurate positioning movements. With a parabolic moving arm, the colposcope is very practical and easy to use. In addition, you can regulate the resistance of individual movements with special locking screws.
- **Specification:**
 - Magnification change stepped in 3 stages. Achievable magnification for f=250mm lens 4x,6x,10x,
 - Stereoscopic tilting head with eyepiece cap regulated in height from: 86 - 126 cm,
 - Axial illumination for shadow-free illumination of the field of view,
 - Micrometric focusing,
 - Eyepieces with angle regulation and dioptric correction,
 - 150W halogen lamp (with cooling) with fiber optics built into the tripod,
 - Green filter,
- **The colposcope is equipped with:**
 - **BRESSER spectacle camera:**
 - 2 MP Full HD,
 - Resolution: 1920 x 1080 px,
 - Ø 23.2 mm cylinder with integrated UV / IR blocking filter,
 - Magnification: approximately 10 x,
 - Image format: jpg, bmp, tif, png,
 - Video format: wmv, avi,
 - **AIO ACL OR-PC 21LP computer:**
 - Windows 10 PL,
 - Video camera software in Polish,
 - Mini IRIS software in trial version (**can be upgraded to full version at additional cost**),
- Mobile, on a 5-arm mobile base with wheels, all with brakes,
- **In the set:**
 - Foot attachment,
 - New, wireless Microsoft keyboard with a touchpad,
 - Power cord,
- Dimensions: 85 x 75 x 135 cm,
- Weight: 48 kg,
- **It has an up-to-date overview and is ready to operate,**
- **Provided with a 12-month valid Technical Passport (Service Report),**
- **Guarantee:**
 - **6 months for the domestic market (Poland),**
 - **3 months for the international market,**
 - **Possible guarantee extension for extra fee, up to 12 months for both, domestic and international market,**
- Financing possibilities (Poland only): Installments, leasing, loan,

